

Alexa Fluor® 647 anti-Cytochrome c Antibody

| | |
|--------------------------|--|
| Catalog# / Size | 612310 / 100 µg |
| Clone | 6H2.B4 |
| Regulatory Status | RUO |
| Other Names | Cyt c |
| Isotype | Mouse IgG1, κ |
| Description | Cytochrome c is a 15 kD protein found in the mitochondrial intermembrane space with a heme-binding domain. Cytochrome c is a component of the electron transport chain; the heme group transfers electrons from cytochrome b-c1 complex to cytochrome oxidase complex. Cytochrome c initiates apoptosis by release to cytoplasm and binding Apaf-1 which activates procaspase 9. Cytochrome c interacts with the cytochrome b-c1 complex, cytochrome oxidase complex, heme, Apaf-1, and Caspase 9 proteins. The 6H2.B4 monoclonal antibody recognizes human, mouse, and rat cytochrome-c and has been shown to be useful for intracellular flow cytometric staining, Western blotting, immunoprecipitation, and immunofluorescence staining. |

Product Details

| | |
|---|---|
| Verified Reactivity | Mouse, Rat, Human |
| Antibody Type | Monoclonal |
| Host Species | Mouse |
| Immunogen | Rat cyt c-OVA |
| Formulation | Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide. |
| Preparation | The antibody was purified by affinity chromatography and conjugated with Alexa Fluor® 647 under optimal conditions. |
| Concentration | 0.5 mg/ml |
| Storage & Handling | The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. Do not freeze. |
| Application | ICFC - Quality tested ICC - Verified |
| Recommended Usage | Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis. For flow cytometric staining, the suggested use of this reagent is ≤0.5 µg per million cells in 100 µl volume. For immunocytochemistry, a concentration range of 2.0 - 5.0 µg/mL is recommended. It is recommended that the reagent be titrated for optimal performance for each application. * Alexa Fluor® 647 has a maximum emission of 668 nm when it is excited at 633 nm / 635 nm. Alexa Fluor® and Pacific Blue™ are trademarks of Life Technologies Corporation. View full statement regarding label licenses |
| Excitation Laser | Red Laser (633 nm) |
| Application Notes | Additional reported applications (for the relevant formats) include: intracellular flow cytometry ⁵ , immunofluorescence microscopy ^{3,5} , immunoprecipitation ⁴ , and immunocytochemistry ⁵ . |
| Application References | 1. Goshorn SC, <i>et al.</i> 1991. <i>J. Biol. Chem.</i> 266:2134. 2. Jemmerson R, <i>et al.</i> 1991. <i>Eur. J. Immunol.</i> 21:143. 3. Chandra D, <i>et al.</i> 2002. <i>J. Biol. Chem.</i> 277:50842. (IF) 4. Semenkova L, <i>et al.</i> 2003. <i>Eur. J. Biochem.</i> 270:4388. (IP) 5. Shih S-F, <i>et al.</i> 2001. <i>J. Biol. Chem.</i> 276:21870. (ICFC ICC IF) 6. She P, <i>et al.</i> 2011. <i>Am J. Physiol Endocrinol Metab.</i> 301:E49. PubMed |
| (PubMed link indicates BioLegend citation) | |

7. McGuire, KA., *et al.* 2011. *J. Virol* 85:10806. [PubMed](#)

Product Citations

1. Sekar G, *et al.* 2022. *iScience*. 25:105064. [PubMed](#)
2. Alcon C, *et al.* 2020. *Cell Death Dis.* 0.8986111111. [PubMed](#)
3. Deng J, *et al.* 2021. *Cell Death Dis.* 12:978. [PubMed](#)
4. Ramsey HE, *et al.* 2018. *Cancer Discov.* 8:1566. [PubMed](#)
5. Meyer LK, *et al.* 2019. *J Clin Invest.* 130:863. [PubMed](#)
6. Olesinski EA, *et al.* 2021. *STAR Protocols.* 2(2):100461. [PubMed](#)
7. Fraser C, *et al.* 2022. *Nat Commun.* 13:5789. [PubMed](#)
8. Jenson J, *et al.* 2017. *Elife.* 10.7554/eLife.25541. [PubMed](#)
9. Bhatt S, *et al.* 2020. *Cancer Cell.* 38(6):872-890.e6. [PubMed](#)
10. Canedo EC, *et al.* 2021. *JCI Insight.* 6: . [PubMed](#)
11. Potter DS, *et al.* 2021. *Cell Death Dis.* 12:741. [PubMed](#)
12. Manzano-Muñoz A, *et al.* 2021. *Front Cell Dev Biol.* 9:695225. [PubMed](#)
13. Alcon C, *et al.* 2021. *Cells.* 10: . [PubMed](#)
14. Victorino F, *et al.* 2021. *eLife.* 0.4166666666666667. [PubMed](#)
15. Dhimolea E, *et al.* 2021. *Cancer Cell.* 39:240. [PubMed](#)
16. Liang Z, *et al.* 2021. *J Immunol.* 207:2039. [PubMed](#)

RRID AB_2565241 (BioLegend Cat. No. 612310)

Antigen Details

| | |
|---------------------------|--|
| Structure | Heme binding domain; 15 kD |
| Distribution | Mitochondrial intermembrane space |
| Function | Component of electron transport chain; heme group transfers electrons from cytochrome b-c1 complex to cytochrome oxidase complex. Initiates apoptosis by release to cytoplasm and binding Apaf-1 which activates procaspase 9 |
| Interaction | Cytochrome b-c1 complex, cytochrome oxidase complex, heme, Apaf-1, Casp9 |
| Biology Area | Apoptosis/Tumor Suppressors/Cell Death, Cell Biology, Mitochondrial Function, Neuroscience, Neuroscience Cell Markers |
| Molecular Family | Mitochondrial Markers |
| Antigen References | <ol style="list-style-type: none">1. Liu X, <i>et al.</i> 1996. <i>Cell.</i> 86:147.2. Li P, <i>et al.</i> 1997. <i>Cell.</i> 91:479.3. Zhang Z, <i>et al.</i> 2003. <i>Gene</i> 312:61.4. Ferguson H, <i>et al.</i> 2003. <i>J. Biol. Chem.</i> 278:45793. |
| Gene ID | 1355 |

Related Protocols


[Surface and Intracellular Cytokine Staining for Flow Cytometry - Video](#)

[Intracellular Flow Cytometry Staining Protocol](#)


Other Formats

Biotin anti-Cytochrome c, FITC anti-Cytochrome c, Purified anti-Cytochrome c, Alexa Fluor® 488 anti-Cytochrome c, Alexa Fluor® 647 anti-Cytochrome c, GMP FITC anti-Cytochrome c

Product Data



C57BL/6 splenocytes were treated with BioLegend's Fixation Buffer and Permeabilization Wash Buffer, and then were stained with Cytochrome C (clone 6H2-B4) Alexa Fluor® 647 (filled histogram) or mouse IgG1, κ Alexa Fluor® 647 isotype control (open histogram).



A-431 cells were fixed with 4% paraformaldehyde for 10 minutes, permeabilized with Triton X-100 for 10 minutes, and blocked with 5% FBS for 60 minutes. Cells were then intracellularly stained with a 1:250 dilution (2.0 µg/mL) of Alexa Fluor® 647 anti-Cytochrome c Antibody, clone 6H2.B4. Nuclei were counterstained with DAPI, and the image was captured with a 60X objective.

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587