

## Alexa Fluor® 700 anti-human HLA-A,B,C Antibody

<b>Catalog# / Size</b>	311437 / 25 tests 311438 / 100 tests
<b>Clone</b>	W6/32
<b>Regulatory Status</b>	RUO
<b>Other Names</b>	Major Histocompatibility Class I, MHC class I
<b>Isotype</b>	Mouse IgG2a, κ
<b>Description</b>	MHC class I antigens associated with β2-microglobulin are expressed by all human nucleated cells. MHC class I molecules are involved in presentation of antigens to CD8 <sup>+</sup> T cells. They play an important role in cell-mediated immune responses and tumor surveillance.

### Product Details

---

<b>Verified Reactivity</b>	Human, Cynomolgus, Rhesus
<b>Reported Reactivity</b>	African Green, Baboon, Cat, Cow, Chimpanzee
<b>Antibody Type</b>	Monoclonal
<b>Host Species</b>	Mouse
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and BSA (origin USA)
<b>Preparation</b>	The antibody was purified by affinity chromatography and conjugated with Alexa Fluor® 700 under optimal conditions.
<b>Concentration</b>	Lot-specific (to obtain lot-specific concentration, please enter the lot number in our <a href="#">Concentration and Expiration Lookup</a> or <a href="#">Certificate of Analysis</a> online tools.)
<b>Storage &amp; Handling</b>	The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. <b>Do not freeze.</b>
<b>Application</b>	<a href="#">FC - Quality tested</a>
<b>Recommended Usage</b>	Each lot of this antibody is quality control tested by <a href="#">immunofluorescent staining with flow cytometric analysis</a> . For flow cytometric staining, the suggested use of this reagent is 5 µl per million cells in 100 µl staining volume or 5 µl per 100 µl of whole blood.  * Alexa Fluor® 700 has a maximum emission of 719 nm when it is excited at 633 nm / 635 nm. Prior to using Alexa Fluor® 700 conjugate for flow cytometric analysis, please verify your flow cytometer's capability of exciting and detecting the fluorochrome.  Alexa Fluor® and Pacific Blue™ are trademarks of Life Technologies Corporation.  <a href="#">View full statement regarding label licenses</a>
<b>Excitation Laser</b>	Red Laser (633 nm)
<b>Application Notes</b>	Clone W6/32 recognizes residues in the N terminus of the human β2-microglobulin molecule <sup>21</sup> .

Additional reported applications (for the relevant formats) include: immunoprecipitation<sup>2</sup>, Western blotting (non-reducing)<sup>3</sup>, immunohistochemical staining of acetone-fixed frozen tissue sections<sup>4,5</sup>, blocking<sup>6,7</sup>, inhibition of NK cell-mediated lysis<sup>10</sup>, and activation<sup>8,9</sup>. Clone W6/32 has been reported not to be suitable for immunohistochemistry on paraffin sections<sup>17</sup>. The LEAF™ purified antibody (Endotoxin < 0.1 EU/µg, Azide-Free, 0.2 µm filtered) is recommended for functional assays. For highly sensitive assays, we recommend Ultra-LEAF™ purified antibody (Cat. No. 311428) with a lower endotoxin limit than standard LEAF™ purified antibodies (Endotoxin < 0.01 EU/µg).

### Application References

(PubMed link indicates BioLegend citation)

1. Darrow TL, *et al.* 1989. *J. Immunol.* 142:3329.
2. Stern P, *et al.* 1987. *J. Immunol.* 138:1088.
3. Tran TM, *et al.* 2001. *Immunogenetics* 53:440.

4. Barbatis C, *et al.* 1981. *Gut* 22:985.
5. Ayyoub M, *et al.* 2004. *Cancer Immunity* 4:7.
6. DeFelice M, *et al.* 1990. *Cell. Immunol.* 126:420.
7. Fayen J, *et al.* 1998. *Int. Immunol.* 10:1347.
8. Turco MC, *et al.* 1988. *J. Immunol.* 141:2275.
9. Geppert TD, *et al.* 1989. *J. Immunol.* 142:3763.
10. Wooden SL, *et al.* 2005. *J. Immunol.* 175:1383.
11. Nagano M, *et al.* 2007. *Blood* 110:151.
12. McLoughlin RM, *et al.* 2008. *J. Immunol.* 181:1323. [PubMed](#)
13. Takahara M, *et al.* 2008. *J. Leukoc. Biol.* 83:742. [PubMed](#)
14. Lunemann A, *et al.* 2008. *J. Immunol.* 181:6170. [PubMed](#)
15. Laing BJ, *et al.* 2010. *J. Thorac Cardiovasc Surg.* 139:1402. [PubMed](#)
16. Yoshino N, *et al.* 2000. *Exp. Anim. (Tokyo)* 49:97. (FC)
17. Vambutas A, *et al.* 2000. *Clin. Diagn. Lab. Immun.* 7:79.
18. Coppieters KT, *et al.* 2012. *J. Exp. Med.* 209:51. (epitope)
19. Crivello P, *et al.* 2013. *Hum Immunol.* 22:100. [PubMed](#)
20. Jung Y, *et al.* 2015. *Mol Cancer Res.* 13:197. [PubMed](#)
21. Shields MJ, Ribaud RK. 1998. *Tissue Antigens.* 51(5):567-70. (epitope)

#### Product Citations

1. Koutsakos M, *et al.* 2019. *Front Immunol.* 10:1158. [PubMed](#)
2. McQueen BE, *et al.* 2020. *Infect Immun.* :88. [PubMed](#)
3. Lim SY, *et al.* 2021. *Biomedicines.* 9:. [PubMed](#)

#### RRID

AB\_2566305 (BioLegend Cat. No. 311437)  
 AB\_2566306 (BioLegend Cat. No. 311438)

### Antigen Details

---

<b>Structure</b>	Ig superfamily
<b>Distribution</b>	All nucleated cells
<b>Function</b>	Antigen presentation
<b>Ligand/Receptor</b>	CD3/TCR, CD8
<b>Biology Area</b>	Immunology, Innate Immunity
<b>Molecular Family</b>	MHC Antigens
<b>Antigen References</b>	1. Barclay AN, <i>et al.</i> Eds. 1993. <i>The Leukocyte Antigen FactsBook.</i> Academic Press Inc. San Diego.
<b>Gene ID</b>	<a href="#">3105</a>

### Related Protocols

---

[Cell Surface Flow Cytometry Staining Protocol](#)


### Other Formats

---

APC anti-human HLA-A,B,C, FITC anti-human HLA-A,B,C, PE anti-human HLA-A,B,C, PE/Cyanine5 anti-human HLA-A,B,C, Purified anti-human HLA-A,B,C, Alexa Fluor® 488 anti-human HLA-A,B,C, Alexa Fluor® 647 anti-human HLA-A,B,C, Pacific Blue™ anti-human HLA-A,B,C, PerCP anti-human HLA-A,B,C, APC/Cyanine7 anti-human HLA-A,B,C, PerCP/Cyanine5.5 anti-human HLA-A,B,C, Ultra-LEAF™ Purified anti-human HLA-A,B,C, PE/Cyanine7 anti-human HLA-A,B,C, Brilliant Violet 510™ anti-human HLA-A,B,C, Alexa Fluor® 700 anti-human HLA-A,B,C, PE/Dazzle™ 594 anti-human HLA-A,B,C, Biotin anti-human HLA-A,B,C, Brilliant Violet 605™ anti-human HLA-A,B,C, APC/Fire™ 750 anti-human HLA-A,B,C, TotalSeq™-A0058 anti-human HLA-A,B,C, TotalSeq™-C0058 anti-human HLA-A,B,C, TotalSeq™-B0058 anti-human HLA-A,B,C, Spark NIR™ 685 anti-human HLA-A,B,C, TotalSeq™-D0058 anti-human HLA-A,B,C

### Product Data

---



Human peripheral blood lymphocytes were stained with HLA-A,B,C (clone W6/32) Alexa Fluor® 700 (filled histogram) or mouse IgG2a, κ Alexa Fluor® 700 isotype control (open histogram).

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587