

## LEGENDplex™ Human CXCL8 (IL-8) Capture Bead B4, 14X

|                          |   |
|--------------------------|---|
| <b>Catalog# / Size</b>   | 741104 / 100 tests (250 µL)   |
| <b>Regulatory Status</b> | RUO   |
| <b>Other Names</b>       | IL8, Interleukin 8, C-X-C Motif Chemokine Ligand 8, Monocyte-Derived Neutrophil Chemotactic Factor, Monocyte-Derived Neutrophil-Activating Peptide, Granulocyte Chemotactic Protein 1, Chemokine (C-X-C Motif) Ligand 8, Beta Endothelial Cell-Derived Neu  |
| <b>Description</b>       | <p>The LEGENDplex™ COVID-19 Cytokine Storm Panel 1 Mix and Match Capture Beads provides capture beads for 14 targets including IL-6, CCL2 (MCP-1), G-CSF, IFN-α2, CCL5 (RANTES), IL-2, IFN-γ, IL-7, IL-1RA, CXCL8 (IL-8), TNF-α, CXCL10 (IP-10), CCL3 (MIP-1a), and IL-10. These beads are recommended for use with each other and the following in a Mix and Match Assay.</p> <p><b>Note:</b> for serum/plasma samples, CCL5 (RANTES) cannot be mixed and matched with other targets and needs to be measured separately due to high dilution requirement.</p> <p>LEGENDplex™ Cytokine Storm Panel 1<br/>Capture Beads are intended for use with the following reagents:<br/>741095 (LEGENDplex™ COVID-19 Cytokine Storm Panel 1 &amp; 2 Standard)<br/>741094 (LEGENDplex™ COVID-19 Cytokine Storm Panel 1 &amp; 2 Detection Antibodies)<br/>741110 (LEGENDplex™ Buffer Set A2)<br/>740377 (LEGENDplex™ Filter Plate)<br/>OR<br/>740379 (LEGENDplex™ V-bottom Plate)</p> |

### Product Details

---

|                            |   |
|----------------------------|---|
| <b>Verified Reactivity</b> | Human   |
| <b>Application</b>         | <a href="#">Multiplex</a>   |
|                            | Learn more about LEGENDplex™ at <a href="https://www.biolegend.com/legendplex">biolegend.com/legendplex</a> |
|                            | Download the <a href="#">LEGENDplex™ software here</a> .  |

### Antigen Details

---

|                |                      |
|----------------|----------------------|
| <b>Gene ID</b> | <a href="#">3576</a> |
|----------------|----------------------|

### Related Protocols

---

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](https://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic,

therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587