

Retrieve-All Antigen Unmasking #2: Basic 10X (Previously Covance catalog# SIG-31924)

Catalog# / Size	928401 / 500 mL
Regulatory Status	RUO
Other Names	Retrieve-All Antigen Unmasking Solution
Previously	Signet Catalog# 1924 Covance Catalog# SIG-31924

Product Details

Formulation	Antigen Unmasking Solution.
Storage & Handling	Store at room temperature.
Application	IHC
Recommended Usage	Each lot of this product is quality control tested by immunohistochemistry. The optimal antigen unmasking condition is a function of temperature (90 to 100°C) and incubation time (5-30 minutes) and should be established by each investigator.
Application Notes	*Concentrated buffers (100x) may show minimal precipitation after long term storage. Such precipitation will not affect product performance. Retrieve-All is an antigen un-masking solution available in the following three pH formulas for use in heat induced un-masking procedures: Retrieve-All 1 (Universal pH 8) Retrieve-All 2 (Basic pH 10) Retrieve-All 3 (Acidic pH 4.8)

Antigen Details

Gene ID	NA
----------------	----

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587