

LEGENDplex™ Buffer Set A

Catalog# / Size 740368 / 100 tests

Regulatory Status RUO

Description LEGENDplex™ Buffer Set A includes the following items:

Setup Beads
SA-PE
Assay Buffer
Wash Buffer
Plate Sealers
Matrix B

This Buffer Set does not include plates or the Capture Beads, Standard, and Detection Antibodies needed to complete a LEGENDplex™ Assay.

These additional components can be purchased separately by selecting the Mix & Match tab at www.biolegend.com/legendplex, and then navigating to your desired assay.

LEGENDplex™ Buffer Set A is designed for use with the following LEGENDplex™ Mix and Match Panels:

Human Anti-Virus Response Panel
Human Essential Immune Response Panel
Human Th Cytokine Panel V02
Mouse Proinflammatory Chemokine Panel
Mouse Th Cytokine Panel V03
NHP Chemokine/Cytokine Panel
Human Diabetes Panel
HU Proinflam. Chemokine Panel 1

Product Details

Application [Multiplex](#)

Learn more about LEGENDplex™ at biolegend.com/legendplex

Download the [LEGENDplex™ software here](#).

Product Citations

1. Aguayo–Mazzucato C, *et al.* 2019. Cell Metab. 30:129. [PubMed](#)
2. Schiferle EB, *et al.* 2021. Sci Adv. 7: . [PubMed](#)

Antigen Details

Gene ID NA

Related Protocols

[LEGENDplex™ Laboratory Protocol - Video](#)

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587