

## FITC anti-human/mouse Granzyme B Antibody

<b>Catalog# / Size</b>	515403 / 25 tests
<b>Clone</b>	GB11
<b>Regulatory Status</b>	RUO
<b>Other Names</b>	Granzyme-2, serine protease B, CCP1, Asp-ase, CTLA-1
<b>Isotype</b>	Mouse IgG1, $\kappa$
<b>Description</b>	Granzyme B is a 32 kD serine protease, also known as granzyme-2, serine protease B, CCP1, Asp-ase, and CTLA-1. Granzyme B is abundantly stored in the granules of cytotoxic T lymphocytes and NK cells. Low level of expression has been reported in granulocytes, B cells, and activated dendritic cells. Granzyme B is crucial for rapid induction of cell death and apoptosis through interaction with mannose-6-phosphate receptor.

### Product Details

---

<b>Verified Reactivity</b>	Human, Mouse
<b>Reported Reactivity</b>	Rat
<b>Antibody Type</b>	Monoclonal
<b>Host Species</b>	Mouse
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and BSA (origin USA)
<b>Preparation</b>	The antibody was purified by affinity chromatography, and conjugated with FITC under optimal conditions.
<b>Concentration</b>	Lot-specific (to obtain lot-specific concentration, please enter the lot number in our <a href="#">Concentration and Expiration Lookup</a> or <a href="#">Certificate of Analysis</a> online tools.)
<b>Storage &amp; Handling</b>	The antibody solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. <b>Do not freeze.</b>
<b>Application</b>	<a href="#">ICFC - Quality tested</a>
<b>Recommended Usage</b>	Each lot of this antibody is quality control tested by <a href="#">intracellular immunofluorescent staining with flow cytometric analysis</a> . For flow cytometric staining, the suggested use of this reagent is 5 $\mu$ l per million cells in 100 $\mu$ l staining volume or 5 $\mu$ l per 100 $\mu$ l of whole blood.
<b>Excitation Laser</b>	Blue Laser (488 nm)
<b>Application References</b>	<ol style="list-style-type: none"><li>1. Wever PC, <i>et al.</i> 1998. <i>Immunology</i>. 93:383</li><li>2. Arens R, <i>et al.</i> 2004. <i>J. Exp. Med.</i> 199:1595</li><li>3. Lima M, <i>et al.</i> 2003. <i>Am. J. Pathol.</i> 163:763</li><li>4. Wiede F, <i>et al.</i> 2014. <i>J Autoimmun.</i> 53:105. <a href="#">PubMed</a></li><li>5. Baker GF, <i>et al.</i> 2014. <i>Cancer Res.</i> 74:5079. <a href="#">PubMed</a></li><li>6. Nacer A, <i>et al.</i> 2014. <i>PLoS Pathog.</i> 10:1004528. <a href="#">PubMed</a></li><li>7. Sharma SK, <i>et al.</i> 2015. <i>J Immunol.</i> 194:5529. <a href="#">PubMed</a></li></ol>

### Product Citations

1. Chan LC, *et al.* 2019. *J Clin Invest.* 129:3324. [PubMed](#)
2. Ajith A, *et al.* 2021. *Front Immunol.* 12:687715. [PubMed](#)
3. Wu MJ, *et al.* 2022. *Cancer Discov.* 12:812. [PubMed](#)
4. AbuEid M, *et al.* 2021. *iScience.* 24(6):102653. [PubMed](#)
5. Donnarumma T, *et al.* 2016. *Cell Rep.* 17:1571-1583. [PubMed](#)
6. Baker G, *et al.* 2014. *Cancer Res.* 74:5079. [PubMed](#)
7. Denk D, *et al.* 2022. *Immunity.* 55:2059. [PubMed](#)
8. Wang W, *et al.* 2022. *World J Gastrointest Oncol.* 14:1124. [PubMed](#)
9. van Vloten JP, *et al.* 2022. *J Immunother Cancer.* 10:. [PubMed](#)
10. Leeansyah E, *et al.* 2015. *PLoS Pathog.* 11: 1005072. [PubMed](#)
11. Shen T, *et al.* 2019. *Cell Death Discov.* 5:123. [PubMed](#)

12. Huang X, *et al.* 2019. *Cancer Immunol Res.* 1.388888889. [PubMed](#)
13. Stegelmeier AA, *et al.* 2022. *Biomedicines.* 10:. [PubMed](#)
14. Zhang S, *et al.* 2022. *Cell Death Dis.* 13:844. [PubMed](#)
15. Ron-Harel N, *et al.* 2019. *Cell Rep.* 28:3011. [PubMed](#)
16. Chow MT *et al.* 2019. *Immunity.* 50(6):1498-1512 . [PubMed](#)
17. Xiao C, *et al.* 2022. *iScience.* 25:103934. [PubMed](#)
18. Figueroa-Romero C, *et al.* 2022. *Front Immunol.* 13:773288. [PubMed](#)
19. Peng Z, *et al.* 2021. *STAR Protocols.* 2(2):100595. [PubMed](#)
20. Di Blasi D, *et al.* 2019. *Cell Mol Gastroenterol Hepatol.* 0.510416667. [PubMed](#)
21. LaFleur MW, *et al.* 2019. *Nat Immunol.* 20:1335. [PubMed](#)
22. Rafiq S, *et al.* 2018. *Nat Biotechnol.* 36:847. [PubMed](#)
23. Kramer KJ, *et al.* 2022. *Nat Commun.* 13:3466. [PubMed](#)
24. Mitchell JE, *et al.* 2021. *Cell Reports.* 35(2):108966. [PubMed](#)
25. Sharma S, *et al.* 2015. *J Immunol.* 194:5529. [PubMed](#)
26. Wu W, *et al.* 2020. *Oncol Lett.* 0.91875. [PubMed](#)
27. Zhang C, *et al.* 2020. *Cell Metabolism.* 31(1):148-161.e5.. [PubMed](#)
28. Lu Y, *et al.* 2020. *Cell.* 180(6):1081-1097. [PubMed](#)
29. Jiang P, *et al.* 2021. *Nat Methods.* 18:1181. [PubMed](#)
30. Islas-Vazquez L, *et al.* 2020. *Biology (Basel).* 9:00. [PubMed](#)
31. Kreutmair S, *et al.* 2021. *Immunity.* . [PubMed](#)
32. Hu J, *et al.* 2020. *Biomed Pharmacother.* 127:110229. [PubMed](#)
33. Lucas C, *et al.* 2021. *Nature.* Online ahead of print. [PubMed](#)
34. Adams NM *et al.* 2018. *Immunity.* 48(6):1172-1182 . [PubMed](#)
35. Wu L, *et al.* 2022. *Theranostics.* 12:842. [PubMed](#)
36. Fujigaki J, *et al.* 2015. *PLoS One.* 10: 0132521. [PubMed](#)
37. Bastiaens G, *et al.* 2016. *Am J Trop Med Hyg.* 94: 663 - 673. [PubMed](#)
38. Wiedemann GM, *et al.* 2020. *Cell Rep.* 33:108498. [PubMed](#)
39. Diaz-Salazar C, *et al.* 2020. *Cell Rep.* 32:108186. [PubMed](#)
40. Huang CH, *et al.* 2020. *Adv Ther (Weinh).* 3:. [PubMed](#)
41. Ulrich-Lewis JT, *et al.* 2022. *Front Immunol.* 13:861710. [PubMed](#)
42. Xu S, *et al.* 2021. *Immunity.* . [PubMed](#)
43. Ringel AE, *et al.* 2020. *Cell.* 183(7):1848-1866.e26. [PubMed](#)
44. Walk J, *et al.* 2020. *Front Immunol.* 2.405555556. [PubMed](#)
45. van Montfoort N, *et al.* 2018. *Cell.* 175:1744. [PubMed](#)
46. Sheng W *et al.* 2018. *Cell.* 174(3):549-563 . [PubMed](#)
47. Vardam-Kaur T, *et al.* 2021. *Oncotarget.* 12:2051. [PubMed](#)
48. Bradley T *et al.* 2018. *Cell.* 175(2):387-399 . [PubMed](#)
49. LaFleur MW, *et al.* 2019. *Nat Commun.* 10:1668. [PubMed](#)
50. Luo Y, *et al.* 2021. *Nat Commun.* 12:3913. [PubMed](#)
51. Webb ER, *et al.* 2022. *iScience.* 25:104995. [PubMed](#)
52. Koikawa K, *et al.* 2021. *Cell.* 184(18):4753-4771.e27. [PubMed](#)
53. Prado-Garcia H, *et al.* 2017. *Cancer Immunol Immunother.* 10.1007/s00262-017-1979-x. [PubMed](#)
54. Simula L *et al.* 2018. *Cell reports.* 25(11):3059-3073 . [PubMed](#)
55. Leclerc M, *et al.* 2019. *Nat Commun.* 10:3345. [PubMed](#)
56. Su S *et al.* 2018. *Cell.* 175(2):442-457 . [PubMed](#)
57. Jeng MY, *et al.* 2018. *J Exp Med.* 215:51. [PubMed](#)
58. Peluso MJ, *et al.* 2021. *Cell Rep.* 36:109518. [PubMed](#)
59. Wiede F, *et al.* 2014. *J Autoimmun.* 53:105. [PubMed](#)
60. Kumagai S, *et al.* 2020. *Immunity.* 53(1):187-203.e8. [PubMed](#)
61. Kmiecik M, *et al.* 2011. *J Immunol.* 187:708. [PubMed](#)
62. Zhang J, *et al.* 2018. *Nature.* 553:91. [PubMed](#)
63. Yu M, *et al.* 2021. *Molecular Cell.* 81(6):1216-1230.e9. [PubMed](#)
64. Maleki KT, *et al.* 2021. *Cell Reports Medicine.* 2(3):100220. [PubMed](#)
65. Zhang H, *et al.* 2021. *Cell Reports.* 35(6):109096. [PubMed](#)
66. Sagebiel AF, *et al.* 2019. *Nat Commun.* 10:975. [PubMed](#)
67. Lettau M, *et al.* 2018. *Int Immunol.* 30:215. [PubMed](#)
68. Oliveira AC *et al.* 2017. *eLife.* 6 pii: e30883. [PubMed](#)
69. Thieme CJ, *et al.* 2020. *Cell Rep Med.* 1:100092. [PubMed](#)
70. Dhar P, *et al.* 2021. *Commun Biol.* 4:905. [PubMed](#)
71. Chen RH, *et al.* 2021. *Front Immunol.* 11:594775. [PubMed](#)
72. Kansy BA, *et al.* 2017. *Cancer Res.* 77:6353. [PubMed](#)
73. Sheppard S, *et al.* 2021. *Cell Reports.* 35(9):109210. [PubMed](#)
74. Hu X, *et al.* 2020. *Neoplasia.* 1.290972222. [PubMed](#)
75. Loo Yau H, *et al.* 2021. *STAR Protocols.* 2(2):100549. [PubMed](#)
76. Loo Yau H, *et al.* 2021. *Molecular Cell.* 81(7):1469-1483.e8. [PubMed](#)
77. Dubrot J, *et al.* 2021. *Immunity.* 54(3):571-585.e6. [PubMed](#)
78. Ajith A, *et al.* 2019. *FASEB J.* 33:5220. [PubMed](#)
79. Geary CD *et al.* 2018. *Cell reports.* 24(8):1949-1957 . [PubMed](#)
80. Lin J, *et al.* 2021. *Cancer Immunol Res.* 9:1413. [PubMed](#)
81. DiMarco AV, *et al.* 2022. *Cancer Immunol Res.* 10:70. [PubMed](#)
82. McAusland TM, *et al.* 2021. *Mol Ther Oncolytics.* 20:306. [PubMed](#)
83. Li Q, *et al.* 2022. *J Biol Chem.* :101753. [PubMed](#)
84. Malenica I, *et al.* 2021. *Nat Commun.* 12:5209. [PubMed](#)
85. Vonderhaar EP, *et al.* 2021. *Cell Mol Gastroenterol Hepatol.* 12:41. [PubMed](#)
86. Ruth JH, *et al.* 2021. *JCI Insight.* 6:e145662. [PubMed](#)
87. Durand J, *et al.* 2015. *J Immunol.* 195: 5035 - 5044. [PubMed](#)
88. Wu L, *et al.* 2020. *Cancer Immunol Res.* 710:8. [PubMed](#)
89. Guo Y, *et al.* 2021. *Nat Immunol.* 22:746. [PubMed](#)
90. Gorby C, *et al.* 2020. *Sci Signal.* :13. [PubMed](#)
91. Aurélien Trompette *et al.* 2018. *Immunity.* 48(5):992-1005 . [PubMed](#)

92. Sethuraman SN, *et al.* 2020. *Theranostics*. 10:3397. [PubMed](#)  
93. Carrion B, *et al.* 2021. *Neurol Neuroimmunol Neuroinflamm*. 8:. [PubMed](#)  
94. Liu Y, *et al.* 2021. *Nat Commun*. 12:6831. [PubMed](#)  
95. Li X, *et al.* 2021. *Front Cell Dev Biol*. 9:647713. [PubMed](#)

**RRID** AB\_2114575 (BioLegend Cat. No. 515403)

## Antigen Details

---

<b>Structure</b>	32 kD serine protease
<b>Distribution</b>	Cytotoxic T-cells and NK cells, low on granulocytes, B cells and activated dendritic cells
<b>Function</b>	Induction of cell death and apoptosis
<b>Ligand/Receptor</b>	Mannose-6-phosphate receptor
<b>Cell Type</b>	B cells, Dendritic cells, NK cells, T cells
<b>Biology Area</b>	Cell Biology, Immunology, Innate Immunity, Neuroscience
<b>Molecular Family</b>	Enzymes and Regulators, Proteases
<b>Antigen References</b>	1. Estebanez-Perpina E, <i>et al.</i> 2000. <i>Biol Chem</i> . 381:1203 2. Griffiths GM. And S. Isaza, <i>et al.</i> 1993. <i>J. Cell Biol</i> . 120:885 3. Spaeny-Dekking EH, <i>et al.</i> 1998. <i>J. Immunol</i> . 160:3610 4. Wagner C, <i>et al.</i> 2008. <i>Mol. Immunol</i> . 45:1761
<b>Gene ID</b>	<a href="#">3002</a> <a href="#">14939</a> <a href="#">171528</a>

## Related Protocols

---

[Surface and Intracellular Cytokine Staining for Flow Cytometry - Video](#)

[Intracellular Flow Cytometry Staining Protocol](#)


## Other Formats

---


FITC anti-human/mouse Granzyme B, Alexa Fluor® 647 anti-human/mouse Granzyme B, Pacific Blue™ anti-human/mouse Granzyme B

## Product Data

---



Human peripheral blood lymphocytes surface stained with CD8 PE, then intracellular stained with GB11 FITC



Human peripheral blood lymphocytes surface stained with CD8 PE, then intracellular stained with mouse IgG1 FITC isotype control

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587