

## APC/Cyanine7 Streptavidin

<b>Catalog# / Size</b>	405208 / 100 µg
<b>Regulatory Status</b>	RUO
<b>Description</b>	Streptavidin binds to biotin with high affinity. Streptavidin-APC/Cyanine7 is useful for detecting biotinylated antibodies. The excitation of APC/Cyanine7 by 600-635 nm laser light induces a fluorescence emission maximum of 780 nm.

### Product Details

---

<b>Verified Reactivity</b>	Human, Mouse, Rat, All Species
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide.
<b>Preparation</b>	Streptavidin is conjugated with APC/Cyanine7 under optimal conditions.
<b>Concentration</b>	0.2 mg/ml (concentration relates to the Streptavidin only component of the conjugate)
<b>Storage &amp; Handling</b>	The Streptavidin APC/Cyanine7 solution should be stored undiluted between 2°C and 8°C, and protected from prolonged exposure to light. <b>Do not freeze.</b>
<b>Application</b>	<a href="#">FC - Quality tested</a> <a href="#">ICFC - Verified</a>
<b>Recommended Usage</b>	Each lot of this Streptavidin APC/Cyanine7 is quality control tested by immunofluorescent staining with flow cytometric analysis. The concentration provided is based upon molecular mass of streptavidin independent of any additional molecular mass that might be added by the APC/Cyanine7 conjugation. For immunofluorescent staining, the recommended use of this reagent is =0.06 µg per 10 <sup>6</sup> cells in 100 µl staining volume. It is recommended that the reagent be titrated for optimal performance for each application.
<b>Excitation Laser</b>	Red Laser (633 nm)
<b>Application Notes</b>	Streptavidin-Allophycocyanin/Cyanine7 (APC/Cyanine7) is useful as a second step reagent for indirect immunofluorescent staining, when used in conjunction with biotinylated primary antibodies. The average molecular weight of Streptavidin-APC/Cyanine7 is 350 kD and Streptavidin alone is 52 kD.
<b>Additional Product Notes</b>	BioLegend is in the process of converting the name APC/Cy7 to APC/Cyanine7. The dye molecule remains the same, so you should expect the same quality and performance from our APC/Cyanine7 products. Please contact <a href="#">Technical Service</a> if you have any questions.

### Application References

(PubMed link indicates BioLegend citation)

- Maresz K, *et al.* 2005. *J. Neurochemistry* 95:437.
- Ponomarev ED, *et al.* 2006. *J. Immunol.* 176:1402. [PubMed](#).
- Ponomarev ED, *et al.* 2007. *J. Immunol.* 178:39.
- Mann MK, *et al.* 2007. *J. Immunol.* 178:3447. [PubMed](#)
- Ponomarev ED, *et al.* 2005. *J. Immunol Methods.* 300:32. [PubMed](#)
- Cao YA, *et al.* 2008. *Blood* [PubMed](#)
- Garidou L, *et al.* 2009. *J. Virol.* [PubMed](#)
- Dolgachev V, *et al.* 2009. *J. Immunol.* 183:5705. [PubMed](#)
- Sountoulidis A, *et al.* 2012. *PLoS One.* 7:e41460. [PubMed](#).
- Izumi T, *et al.* 2013. *Cytotherapy.* 15:481. [PubMed](#)
- Hamanaka S, *et al.* 2013. *Biochem Biophys Res Commun.* 435:586. [PubMed](#)
- Prince J, *et al.* 2013. *J. Immunol.* 191:4709. [PubMed](#)
- Phanse Y, *et al.* 2013. *Acta Biomater.* 9:8902. [PubMed](#)
- Orre M, *et al.* 2014. *Neurobiol Aging.* 35:1. [PubMed](#)
- Barbera M, *et al.* 2014. *Gut.* 64:11. [PubMed](#)
- Zhang Y, *et al.* 2014. *Infect Immun.* 82:1698. [PubMed](#)
- Peterson BC, *et al.* 2014. *J Leukoc Biol.* 95:809. [PubMed](#)
- Barbera M, *et al.* 2015. *Gut.* 64:11. [PubMed](#)
- Charmsaz S, *et al.* 2015. *PLoS One.* 10:130692. [PubMed](#)

### Product Citations

- Hutter K, *et al.* 2022. *Front Immunol.* 13:967914. [PubMed](#)
- Zhang Y, *et al.* 2014. *Infect Immun.* 82:1689. [PubMed](#)

3. Petersen B, *et al.* 2014. J Leukoc Biol. 95:809. [PubMed](#)
4. Ponomarev E, *et al.* 2006. J Immunol. 176:1402. [PubMed](#)
5. Wang L, *et al.* 2019. Cell Rep. 29:1848. [PubMed](#)
6. Schloss MJ, *et al.* 2022. Nat Immunol. 23:605. [PubMed](#)
7. Rohde D, *et al.* 2022. Nat Cardiovasc Res. 1:28. [PubMed](#)
8. Marco Barros R, *et al.* 2016. Cell. 167: 203-218. [PubMed](#)
9. Gudjonsson A, *et al.* 2019. Front Immunol. 10:1086. [PubMed](#)
10. Gutiérrez-Gutiérrez Ó, *et al.* 2021. EMBO Rep. 22:e52905. [PubMed](#)
11. Orre M, *et al.* 2014. Neurobiol Aging. 35:1. [PubMed](#)
12. Bonnet C, *et al.* 2021. iScience. 24:103399. [PubMed](#)
13. Barbera M, *et al.* 2015. Gut. 64:11. [PubMed](#)
14. Wu T, *et al.* 2021. BMC Biotechnol. 21:42. [PubMed](#)
15. Sánchez-Salguero E, *et al.* 2021. Front Immunol. 12:712130. [PubMed](#)
16. Tu X, *et al.* 2022. Nat Commun. 13:6977. [PubMed](#)
17. Prince J, *et al.* 2013. J Immunol. 191:4709. [PubMed](#)
18. Phanse Y, *et al.* 2013. Acta Biomater. 9:8902. [PubMed](#)
19. Seri Y, *et al.* 2015. Sci Rep. 5: 13041. [PubMed](#)
20. Lin K, *et al.* 2016. Sci Rep. 6:31444. [PubMed](#)
21. Kimball A, *et al.* 2018. Arterioscler Thromb Vasc Biol. 38:1102. [PubMed](#)
22. Johansen AZ, *et al.* 2022. Pharmaceuticals. 14:1. [PubMed](#)
23. Keller EJ, *et al.* 2021. Front Immunol. 12:616064. [PubMed](#)
24. Yeo S, *et al.* 2016. Cancer Res. 76: 3397 - 3410. [PubMed](#)
25. Izumi T, *et al.* 2013. Cytotherapy. 15:481. [PubMed](#)
26. Ponomarev ED, *et al.* 2005. J Immunol Methods. 300:32:00. [PubMed](#)
27. Charmsaz S, *et al.* 2015. PLoS One. 10:130692. [PubMed](#)
28. Ritter M, *et al.* 2020. Ann Hematol. 99:2329. [PubMed](#)
29. Barman PK, *et al.* 2019. J Immunol. 202:2720. [PubMed](#)
30. Greenwald AC, *et al.* 2019. J Exp Med. 216:215. [PubMed](#)
31. Zhu H, *et al.* 2020. J Bone Miner Res. 35:698. [PubMed](#)
32. Chen Y, *et al.* 2015. PLoS One. 10: 0135217. [PubMed](#)
33. Frodermann V, *et al.* 2019. Nat Med. 25:1761. [PubMed](#)
34. Harnett MM, *et al.* 2022. Front Immunol. 13:953053. [PubMed](#)
35. Omer OS, *et al.* 2020. Methods Mol Biol. 2121:199. [PubMed](#)
36. Lebert-Ghali C, *et al.* 2016. Blood. 127: 87-90. [PubMed](#)
37. Dolgachev V, *et al.* 2009. J Immunol. 183:5705. [PubMed](#)
38. Delacher M, *et al.* 2020. Immunity. 52:295. [PubMed](#)
39. Chen Z, *et al.* 2019. J Exp Med. 216:152. [PubMed](#)
40. Freter R, *et al.* 2021. Nat Commun. 12:3318. [PubMed](#)
41. Shimizu T, *et al.* 2016. J Exp Med. 213: 1479 - 1496. [PubMed](#)
42. Gillespie A, *et al.* 2016. PLoS One. 12: 1005419. [PubMed](#)
43. Chao D, *et al.* 2015. PLoS One. 10: 0138107. [PubMed](#)
44. Hamanaka S, *et al.* 2013. Biochem Biophys Res Commun. 435:586. [PubMed](#)
45. Mann M, *et al.* 2007. J Immunol. 178:3447. [PubMed](#)
46. Cao Y, *et al.* 2008. Blood. 112:4494. [PubMed](#)
47. Kulkarni N, *et al.* 2017. Eur J Appl Physiol. 10.1007/s00421-017-3704-z. [PubMed](#)
48. Keller EJ, *et al.* 2022. Int J Mol Sci. 23:1. [PubMed](#)
49. King R, *et al.* 2021. Sci Adv. 7:eabj5293. [PubMed](#)
50. Ho JY, *et al.* 2021. Mol Ther Methods Clin Dev. 21:237. [PubMed](#)
51. Girardi R, *et al.* 2015. Nat Commun. 6: 8487. [PubMed](#)
52. Chen Z, *et al.* 2020. Cell Death Dis. 11:180. [PubMed](#)
53. Sountoulidis A, *et al.* 2012. PLoS One. 7:e41460. [PubMed](#)
54. Kimball AS, *et al.* 2020. Thromb Haemost. 120:289. [PubMed](#)
55. Thürmer M, *et al.* 2022. Nat Commun. 13:2982. [PubMed](#)
56. Mesnieres M, *et al.* 2021. Cell Reports. 36(8):109618. [PubMed](#)
57. Wang J, *et al.* 2022. Elife. 11:1. [PubMed](#)
58. Luo H, *et al.* 2019. Cell Rep. 26:945. [PubMed](#)
59. Garidou L, *et al.* 2009. J Virol. 83:8905. [PubMed](#)
60. Schafer JB, *et al.* 2022. J Biol Chem. 297:102694. [PubMed](#)
61. Furusawa J, *et al.* 2016. PLoS Pathog. 12: 1005507. [PubMed](#)
62. Sen D, *et al.* 2016. PLoS One. 11:e0165064. [PubMed](#)
63. Li ZG, *et al.* 2020. Chin Med J (Engl). 1943:133. [PubMed](#)
64. Song W, *et al.* 2022. Immunity. 55:290. [PubMed](#)
65. Man N, *et al.* 2021. JCI Insight. 6:1. [PubMed](#)

## Antigen Details

---

Gene ID                      NA

## Related Protocols

---

[Cell Surface Flow Cytometry Staining Protocol](#)

## Intracellular Flow Cytometry Staining Protocol

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587