



GMP Recombinant Human IL-7 (carrier-free)

Catalog# / Size 581914 / 25 μg 581916 / 100 μg

Other Names Interleukin 7, Lymphopoietin-1 (LP-1), IL7

Description IL-7 was initially described as a stromal derived factor which is capable of

inducing the growth of pre-B cells in vitro. IL-7 acts on a variety of cells through its receptor (IL-7R), a heterodimer consisting of IL-7Rα (CD127) and a common γ chain (CD132, yc) shared by other cytokine (IL-2, IL-4, IL-9, IL-15, and IL-21) receptors. In addition, IL-7Rα is shared with TSLP. The generation of IL-7deficient and IL-7Rα-deficient mice and monoclonal antibody blocking experiments confirmed the requirement of IL-7 for B-cell development in mice. Nevertheless, mutations in the α chain of the IL-7 receptor in patients with severe combined immunodeficiency (SCID) confirmed that IL-7 is indispensable for T-cell development in humans. However, the presence of B cells in these individuals suggests important differences between the role of IL-7 in murine and human lymphocyte development. Thus, although human B-cell development does not appear to require IL-7, immature human B cells do proliferate in response to IL-7. Nevertheless, most recent information suggests that IL-7 dependence in human lymphopoiesis increases during the progression of ontogeny in cord blood and bone marrow. IL-7 can be associated to hepatocyte growth factor (HGFβ) to form a hybrid cytokine (IL-7/HGFβ), which induces greater proliferation of CFU-S, SLPs, and pre-pro-B cells than does native IL-7. The hybrid cytokine signals through both IL-7R (IL-7Rα plus γc) and c-Met. IL-7 has a thymic antiapoptotic effect and induces the expression of antiapoptotic proteins Bcl-2 and Mcl-1 and the inhibition of proapoptotic proteins Bax and Bad. In addition, IL-7 is a key regulator of glucose uptake in T lymphocytes. TGF-β has been shown to down-regulate IL-7 mRNA and protein secretion from human bone marrow stromal cells. In addition, TGF- β inhibits IL-7-induced proliferation of pre-B cells.

Product Details

Source Human IL-7, amino acids Asp26-His177 (Accession # NP_000871.1) was

expressed in 293E cells.

Molecular Mass The 165 amino acid recombinant protein has a predicted molecular mass of

approximately 18.8 kD. The DTT-reduced and non-reduced protein migrate at

approximately 20 - 30 kD by SDS-PAGE.

Purity > 95%, as determined by Coomassie stained SDS-PAGE

Formulation 0.1 µm filtered protein solution is in PBS, pH 7.2.

Concentration 500 µg/mL

Storage & Handling Unopened vial can be stored between 2°C and 8°C for up to 2 weeks, at -20°C for

up to six months, or at -70°C or colder until the expiration date. For maximum results, quick spin vial prior to opening. The protein can be aliquoted and stored at -20°C or colder. Stock solutions can also be prepared at 50 - 100 µg/mL in appropriate sterile buffer, carrier protein such as 0.2 - 1% endotoxin-free BSA or HSA can be added when preparing the stock solution. Aliquots can be stored between 2°C and 8°C for up to one week or stored at -20°C or colder for up to 3

months. Avoid repeated freeze/thaw cycles.

Activity ED₅₀ = 0.2 - 1.2 ng/mL, as determined by a dose-dependent stimulation in a

proliferation assay with PHA-activated human peripheral blood lymphocytes (PBL).

Deep Blue Cell Viability™ Kit (Cat. No. 424701) is used to measure the

proliferation.

Application Bioassay

Cell Culture

Application Notes BioLegend carrier-free recombinant proteins provided in liquid format are shipped

on blue ice. Our comparison testing data indicates that when handled and stored as recommended, the liquid format has equal or better stability and shelf-life compared to commercially available lyophilized proteins after reconstitution. Our liquid proteins are verified in-house to maintain activity after shipping on blue ice and are backed by our 100% satisfaction guarantee. If you have any concerns,

contact us at tech@biolegend.com.

Disclaimer

GMP Recombinant Proteins. BioLegend GMP recombinant proteins are manufactured in a dedicated GMP facility and compliant with ISO 13485:2016. For research or ex vivo cell processing use. Not for use in diagnostic or therapeutic procedures. Our processes include:

- · Batch-to-batch consistency
- Material traceability
- Documented procedures
- · Documented employee training
- · Equipment maintenance and monitoring records
- · Lot-specific certificates of analysis
- Quality audits per ISO 13485:2016
- · QA review of releaed products

BioLegend GMP recombinant proteins are manufactured and tested in accordance with USP Chapter 1043, Ancillary Materials for Cell, Gene and Tissue-Engineered Products and Ph. Eur. Chapter 5.2.12.

Antigen Details

Distribution IL-7 is produced by epithelial cells in thymus, bone marrow, and intestine.

Additional sites of IL-7 production include epithelial goblet cells, keratinocytes, fetal liver, adult liver, dendritic cells, skeletal muscle cells, fibroblastic reticular cells, and

follicular dendritic cells.

Function IL-7 induces proliferation of human immature B cells, and it is critical for T-cell

development and peripheral T-cell homeostasis. Regulates T folicular helper cells.

Interaction B cells, Hematopoietic stem and progenitors

Ligand/Receptor IL-7Rα (CD127), common gamma chain (yc or CD132)

Bioactivity Measured by its ability to induce proliferation of 7TD1 cell

Cell Type B cells, Dendritic cells, Epithelial cells, Leukemia, Lymphocytes, NK cells, T cells,

Tfh, Thymocytes


Biology Area Cell Biology, Cell Proliferation and Viability, Immuno-Oncology, Immunology, Stem

Molecular Family Cytokine/Chemokine Receptors, Cytokines/Chemokines

Antigen References

- 1. Park LS, et al. 1990. J. Exp. Med. 171:1073.
- 2. Fry TJ and Mackall CL. 2002. Blood 99:3892.
- Lai L, et al. 2006. Blood 107:1776.
 Link A, et al. 2007. Nat. Immunol. 8:1255.
- 5. Wofford J, et al. 2008. Blood 111:2101.
- 6. Parrish YK, et al. 2009. J. Immunol. 182:4255.
- 7. Saini M, et al. 2009. Blood 113:5793.
- 8. Guimond M, et al. 2009. Nat. Immunol. 10:149.

Gene ID 3574



GMP recombinant human IL-7 induces proliferation of PHA-activated human peripheral blood lymphocytes (PBL) in a dose-dependent manner with an ED₅₀ range of 0.2 - 1.2 ng/mL.

Symbols Glossary*

Symbol	Meaning	Symbol Title	Symbol No.	Symbol	Meaning	Symbol Title	Symbol No.
REF	Catalog number	Catalogue number	5.1.6	(i	Indicates the need for the user to consult the instructions for use.	Consult instructions for use	5.4.3
1	Indicates the temperature limits to which the medical device can be safely exposed.	Temperature limit	5.3.7	类	Indicates a medical device that needs protection from light sources.	Keep away from sunlight	5.3.2
X	Indicates the upper limit of temperature to which the medical device can be safely exposed.	temperature	5.3.6	Ω	Indicates the date after which the medical device is not to be used.	Use-by date	5.1.4
4	Indicates the medical device manufacturer.	Manufacturer	5.1.1	EC REP	Indicates the authorized representative in the European Community.	Authorized representative in the European Community	
LOT	Indicates the manufacturer's batch code so that the batch or lot can be identified.	Batch code	5.1.5	IVD	Indicates a medical device that is intended to be used as an in vitro diagnostic medical device.	In vitro diagnostic medical device	5.5.1

^{*} Symbol information is from EN ISO 15223-1:2016 Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements

For research and *ex vivo* cell processing use. Not for diagnostic or therapeutic procedures. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/terms). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to reverse engineer functionally similar materials without written approval of BioLegend. By use of these products you accept the terms and conditions of any and all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research and ex vivo cell processing use only and not intended for human or animal diagnostic or therapeutic use.

8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587