

## TotalSeq™-C0051 anti-human CD14 Antibody

<b>Catalog# / Size</b>	367137 / 10 µg
<b>Clone</b>	63D3
<b>Regulatory Status</b>	RUO
<b>Other Names</b>	Monocyte differentiation antigen CD14, myeloid cell-specific leucine-rich glycoprotein, LPS receptor
<b>Isotype</b>	Mouse IgG1, κ
<b>Barcode Sequence</b>	CAATCAGACCTATGA
<b>Description</b>	CD14 is a 53-55 kD glycosylphosphatidylinositol (GPI)-linked membrane glycoprotein that is also known as the LPS receptor. CD14 is expressed at high levels on monocytes and macrophages, and at lower levels on granulocytes. Some dendritic cell populations such as interfollicular dendritic cells, reticular dendritic cells, and Langerhans cells have also been reported to express CD14. As a high-affinity receptor for LPS, CD14 is involved in the clearance of gram-negative pathogens and in the upregulation of adhesion molecules and cytokine expression in monocytes and neutrophils.

### Product Details

---

<b>Verified Reactivity</b>	Human
<b>Antibody Type</b>	Monoclonal
<b>Host Species</b>	Mouse
<b>Immunogen</b>	Purified human peripheral blood monocytes.
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 1 mM EDTA.
<b>Preparation</b>	The antibody was purified by chromatography and conjugated with TotalSeq™-C oligomer under optimal conditions.
<b>Concentration</b>	0.5 mg/ml
<b>Storage &amp; Handling</b>	The antibody solution should be stored undiluted between 2°C and 8°C. <b>Do not freeze.</b>
<b>Application</b>	<a href="#">PG - Quality tested</a>
<b>Recommended Usage</b>	<p>Each lot of this antibody is quality control tested by <a href="#">immunofluorescent staining with flow cytometric analysis</a> and the oligomer sequence is confirmed by sequencing. TotalSeq™-C antibodies are compatible with 10x Genomics Chromium Single Cell Immune Profiling <a href="#">Solution</a>.</p> <p>To maximize performance, it is strongly recommended that the reagent be titrated for each application, and that you centrifuge the antibody dilution before adding to the cells at 14,000xg at 2 - 8°C for 10 minutes. Carefully pipette out the liquid avoiding the bottom of the tube and add to the cell suspension. For Proteogenomics analysis, the suggested starting amount of this reagent for titration is ≤ 1.0 µg per million cells in 100 µL volume. Refer to the corresponding TotalSeq™ protocol for specific staining instructions.</p> <p>Buyer is solely responsible for determining whether Buyer has all intellectual property rights that are necessary for Buyer's intended uses of the BioLegend TotalSeq™ products. For example, for any technology platform Buyer uses with TotalSeq™, it is Buyer's sole responsibility to determine whether it has all necessary third party intellectual property rights to use that platform and TotalSeq™ with that platform.</p>
<b>Additional Product Notes</b>	<p>TotalSeq™ reagents are designed to profile protein levels at a single cell level following an optimized protocol similar to the CITE-seq workflow. A compatible single cell device (e.g. <a href="#">10x Genomics Chromium System and Reagents</a>) and sequencer (e.g. Illumina analyzers) are required. Please contact <a href="#">technical support</a> for more information, or visit <a href="http://biolegend.com/totalseq">biolegend.com/totalseq</a>.</p> <p>The barcode flanking sequences are CGGAGATGTGTATAAGAGACAGNNNNNNNNNN (PCR handle), and NNNNNNNNNCCCATATAAGA*A*A (capture sequence). N represents either randomly</p>

selected A, C, G, or T, and \* indicates a phosphorothioated bond, to prevent nuclease degradation.

View more applications data for this product in our [Scientific Poster Library](#).

## Application References

(PubMed link indicates BioLegend citation)

1. Fridlender ZG, *et al.* 1999. *Hum. Immunol.* 11:1028.
2. Devitt A, *et al.* 1998. *Nature* 6675:505.

RRID AB\_2810576 (BioLegend Cat. No. 367137)

## Antigen Details

---

<b>Structure</b>	GPI-linked membrane glycoprotein.
<b>Distribution</b>	Monocytes, macrophages, dendritic cells, and granulocytes.
<b>Function</b>	LPS receptor, clearance of Gram-negative pathogens.
<b>Interaction</b>	LPS.
<b>Ligand/Receptor</b>	LPS.
<b>Cell Type</b>	Dendritic cells, Granulocytes, Macrophages, Monocytes, Neutrophils
<b>Biology Area</b>	Cell Biology, Immunology, Neuroinflammation, Neuroscience
<b>Molecular Family</b>	Adhesion Molecules, CD Molecules

## Antigen References

1. Stocks SC, *et al.* 1990. *Biochem. J.* 268:275.
2. Wright SD, *et al.* 1990. *Science* 4975:1431.

Gene ID [929](#)

## Related Protocols

---

[TotalSeq™-B or -C with 10x Feature Barcoding Technology](#)

## Other Formats

---

Purified anti-human CD14, PE anti-human CD14, Biotin anti-human CD14, APC/Cyanine7 anti-human CD14, PE/Cyanine7 anti-human CD14, PerCP/Cyanine5.5 anti-human CD14, Alexa Fluor® 700 anti-human CD14, APC anti-human CD14, FITC anti-human CD14, APC/Fire™ 750 anti-human CD14, Pacific Blue™ anti-human CD14, Brilliant Violet 510™ anti-human CD14, Brilliant Violet 605™ anti-human CD14, Alexa Fluor® 647 anti-human CD14, Alexa Fluor® 488 anti-human CD14, TotalSeq™-A0051 anti-human CD14, PE/Dazzle™ 594 anti-human CD14, Brilliant Violet 750™ anti-human CD14, TotalSeq™-C0051 anti-human CD14, Brilliant Violet 421™ anti-human CD14, Brilliant Violet 785™ anti-human CD14, Brilliant Violet 711™ anti-human CD14, TotalSeq™-B0051 anti-human CD14, Spark Blue™ 550 anti-human CD14, Spark NIR™ 685 anti-human CD14, PE/Fire™ 640 anti-human CD14, APC/Fire™ 810 anti-human CD14, PerCP anti-human CD14, PE/Fire™ 700 anti-human CD14

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587

