

TotalSeq™-A0064 anti-human CD123 Antibody

Catalog# / Size	306037 / 10 µg
Clone	6H6
Regulatory Status	RUO
Other Names	IL-3R α , IL-3 Receptor alpha
Isotype	Mouse IgG1, κ
Barcode Sequence	CTTCACTCTGTCAGG
Description	CD123 is the 70 kD transmembrane α chain of the IL-3 receptor. Alone, CD123 binds IL-3 with low affinity; when CD123 associates with CDw131 (common β chain), it binds IL-3 with high affinity. CD123 does not transduce intracellular signals upon binding IL-3 and requires the β chain for this function. CD123 is expressed by myeloid precursors, macrophages, dendritic cells, mast cells, basophils, megakaryocytes, and some B cells.

Product Details

Verified Reactivity	Human
Reported Reactivity	Rhesus
Antibody Type	Monoclonal
Host Species	Mouse
Immunogen	Human IL-3R α transfected COS cells.
Formulation	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 1 mM EDTA.
Preparation	The antibody was purified by chromatography and conjugated with TotalSeq™-A oligomer under optimal conditions.
Concentration	0.5 mg/ml
Storage & Handling	The antibody solution should be stored undiluted between 2°C and 8°C. Do not freeze.
Application	PG - Quality tested
Recommended Usage	<p>Each lot of this antibody is quality control tested by immunofluorescent staining with flow cytometric analysis and the oligomer sequence is confirmed by sequencing. TotalSeq™-A antibodies are compatible with 10x Genomics Single Cell Gene Expression Solutions.</p> <p>To maximize performance, it is strongly recommended that the reagent be titrated for each application, and that you centrifuge the antibody dilution before adding to the cells at 14,000xg at 2 - 8°C for 10 minutes. Carefully pipette out the liquid avoiding the bottom of the tube and add to the cell suspension. For Proteogenomics analysis, the suggested starting amount of this reagent for titration is ≤ 1.0 µg per million cells in 100 µL volume. Refer to the corresponding TotalSeq™ protocol for specific staining instructions.</p> <p>Buyer is solely responsible for determining whether Buyer has all intellectual property rights that are necessary for Buyer's intended uses of the BioLegend TotalSeq™ products. For example, for any technology platform Buyer uses with TotalSeq™, it is Buyer's sole responsibility to determine whether it has all necessary third party intellectual property rights to use that platform and TotalSeq™ with that platform.</p>
Application Notes	Clone 6H6 does not inhibit IL-3 binding to low- or high-affinity IL-3Rs. Additional reported applications (for the relevant formats) include: Western blotting ¹ , immunoprecipitation ¹ , and immunohistochemical staining of acetone-fixed frozen sections ² and also paraformaldehyde fixed paraffin embedded tissue ⁷ .
Additional Product Notes	TotalSeq™ reagents are designed to profile protein levels at a single cell level following an optimized protocol similar to the CITE-seq workflow. A compatible single cell device (e.g. 10x Genomics Chromium System and Reagents) and sequencer (e.g. Illumina analyzers) are required.

Please contact [technical support](#) for more information, or visit biolegend.com/totalseq.

The barcode flanking sequences are CCTTGGCACCCGAGAATTCCA (PCR handle), and BAAA*A*A (capture sequence). B represents either C, G, or T, and * indicates a phosphorothioated bond, to prevent nuclease degradation.

View more applications data for this product in our [Scientific Poster Library](#).

Application References

(PubMed link indicates BioLegend citation)

1. Sun Q, *et al.* 1996. *Blood* 87:83. (IP, WB)
2. Herling M, *et al.* 2003. *Blood* 101:5007. (IHC)
3. Charles N, *et al.* 2010. *Nat. Med.* 16:701. (FC) [PubMed](#)
4. Martin-Gayo E, *et al.* 2010. *Blood* 115:5366. [PubMed](#)
5. Chen SC, *et al.* 2010. *Arch Dermatol Res.* 302:113. [PubMed](#)
6. Liu Y, *et al.* 2012. *Food Chem Toxicol.* 50:1920. [PubMed](#)
7. Peduzzi E, *et al.* 2007. *J. Invest. Dermatol.* 127:638. (IHC)

Product Citations

1. Cook CP, *et al.* 2022. *Cell Rep Med.* 3:100715. [PubMed](#)
2. Stuart T, *et al.* 2019. *Cell.* 177:1888. [PubMed](#)
3. Witkowski MT, *et al.* 2020. *Cancer Cell.* 37:867. [PubMed](#)
4. Swanson E, *et al.* 2021. *eLife.* 10:00. [PubMed](#)
5. Hao Y, *et al.* 2021. *Cell.* 184:3573. [PubMed](#)

RRID

AB_2749977 (BioLegend Cat. No. 306037)

Antigen Details

Structure	Ig superfamily, type I transmembrane glycoprotein, associates with CDw131, 70 kD
Distribution	Myeloid precursors, basophils, mast cells, macrophages, dendritic cells, megakaryocytes, subset of lymphocytes
Function	Hematopoietic cell proliferation, differentiation
Ligand/Receptor	IL-3
Cell Type	Basophils, Dendritic cells, Hematopoietic stem and progenitors, Lymphocytes, Macrophages, Mast cells, Megakaryocytes
Biology Area	Immunology
Molecular Family	CD Molecules, Cytokine/Chemokine Receptors
Antigen References	1. Miyajima A, <i>et al.</i> 1993. <i>Blood</i> 82:1960.
Gene ID	3563

Related Protocols

[TotalSeq™-A Antibodies and Cell Hashing with 10x Single Cell 3' Reagent Kit v3 3.1 Protocol](#)

Other Formats

Biotin anti-human CD123, PE anti-human CD123, Purified anti-human CD123, PE/Cyanine5 anti-human CD123, PE/Cyanine7 anti-human CD123, APC anti-human CD123, FITC anti-human CD123, PerCP/Cyanine5.5 anti-human CD123, Brilliant Violet 421™ anti-human CD123, Brilliant Violet 650™ anti-human CD123, Brilliant Violet 510™ anti-human CD123, Alexa Fluor® 647 anti-human CD123, Brilliant Violet 605™ anti-human CD123, Purified anti-human CD123 (Maxpar® Ready), Brilliant Violet 711™ anti-human CD123, Brilliant Violet 785™ anti-human CD123, PE/Dazzle™ 594 anti-human CD123, Alexa Fluor® 488 anti-human CD123, TotalSeq™-A0064 anti-human CD123, Alexa Fluor® 700 anti-human CD123, APC/Fire™ 750 anti-human CD123, Pacific Blue™ anti-human CD123, TotalSeq™-C0064 anti-human CD123, TotalSeq™-B0064 anti-human CD123, TotalSeq™-D0064 anti-human CD123, PerCP anti-human CD123, GMP PE/Cyanine7 anti-human CD123, APC/Fire™ 810 anti-human CD123 Antibody

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587