

## TotalSeq™-A0392 anti-human CD15 (SSEA-1) Antibody

<b>Catalog# / Size</b>	323046 / 10 µg
<b>Clone</b>	W6D3
<b>Regulatory Status</b>	RUO
<b>Workshop</b>	HCDM listed
<b>Other Names</b>	Lewis X, 3-FAL, 3-FL, LNFP III, Lex, SSEA-1, X-hapten
<b>Isotype</b>	Mouse IgG1, κ
<b>Barcode Sequence</b>	TCACCAAGTACCTAGT
<b>Description</b>	CD15 is 3-fucosyl-N-acetyllactosamine (3-FAL) also known as Lewis X, 3-FAL, X-hapten, and SSEA-1. CD15 is expressed on granulocytes and monocytes. It has also been shown to be expressed on Langerhans cells. CD15 has been implicated in adhesion as well as chemotaxis, phagocytosis, and bactericidal activity.

### Product Details

<b>Verified Reactivity</b>	Human
<b>Antibody Type</b>	Monoclonal
<b>Host Species</b>	Mouse
<b>Immunogen</b>	WERI-RB-1 retinoblastoma cell line
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 1 mM EDTA.
<b>Preparation</b>	The antibody was purified by chromatography and conjugated with TotalSeq™-A oligomer under optimal conditions.
<b>Concentration</b>	0.5 mg/ml
<b>Storage &amp; Handling</b>	The antibody solution should be stored undiluted between 2°C and 8°C. <b>Do not freeze.</b>
<b>Application</b>	<a href="#">PG - Quality tested</a>
<b>Recommended Usage</b>	<p>Each lot of this antibody is quality control tested by <a href="#">immunofluorescent staining with flow cytometric analysis</a> and the oligomer sequence is confirmed by sequencing. TotalSeq™-A antibodies are compatible with 10x Genomics Single Cell Gene Expression <a href="#">Solutions</a>.</p> <p>To maximize performance, it is strongly recommended that the reagent be titrated for each application, and that you centrifuge the antibody dilution before adding to the cells at 14,000xg at 2 - 8°C for 10 minutes. Carefully pipette out the liquid avoiding the bottom of the tube and add to the cell suspension. For Proteogenomics analysis, the suggested starting amount of this reagent for titration is ≤ 1.0 µg per million cells in 100 µL volume. Refer to the corresponding TotalSeq™ protocol for specific staining instructions.</p> <p>Buyer is solely responsible for determining whether Buyer has all intellectual property rights that are necessary for Buyer's intended uses of the BioLegend TotalSeq™ products. For example, for any technology platform Buyer uses with TotalSeq™, it is Buyer's sole responsibility to determine whether it has all necessary third party intellectual property rights to use that platform and TotalSeq™ with that platform.</p>
<b>Additional Product Notes</b>	<p>TotalSeq™ reagents are designed to profile protein levels at a single cell level following an optimized protocol similar to the CITE-seq workflow. A compatible single cell device (e.g. <a href="#">10x Genomics Chromium System and Reagents</a>) and sequencer (e.g. Illumina analyzers) are required. Please contact <a href="#">technical support</a> for more information, or visit <a href="#">biolegend.com/totalseq</a>.</p> <p>The barcode flanking sequences are CCTTGGCACCCGAGAATTCCA (PCR handle), and BAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA*A*A (capture sequence). B represents either C, G, or T, and * indicates a phosphorothioated bond, to prevent nuclease degradation.</p>

View more applications data for this product in our [Scientific Poster Library](#).

## Product Citations

1. Witkowski MT, *et al.* 2020. *Cancer Cell*. 37:867. [PubMed](#)
2. Hao Y, *et al.* 2021. *Cell*. 184:3573. [PubMed](#)

## RRID

AB\_2734304 (BioLegend Cat. No. 323046)

## Antigen Details

---

<b>Structure</b>	Poly-N-acetyllactosamine
<b>Distribution</b>	Neutrophils, eosinophils, monocytes
<b>Function</b>	Adhesion
<b>Cell Type</b>	Embryonic Stem Cells, Eosinophils, Neutrophils
<b>Biology Area</b>	Cell Biology, Immunology, Neuroscience, Neuroscience Cell Markers, Stem Cells
<b>Molecular Family</b>	CD Molecules
<b>Antigen References</b>	1. Stocks SC, <i>et al.</i> 1990. <i>Biochem. J.</i> 268:275.
<b>Gene ID</b>	<a href="#">2526</a>

## Related Protocols

---

[TotalSeq™-A Antibodies and Cell Hashing with 10x Single Cell 3' Reagent Kit v3 3.1 Protocol](#)

## Other Formats

---

Brilliant Violet 785™ anti-human CD15 (SSEA-1), Purified anti-human CD15 (SSEA-1), FITC anti-human CD15 (SSEA-1), PE anti-human CD15 (SSEA-1), APC anti-human CD15 (SSEA-1), Alexa Fluor® 488 anti-human CD15 (SSEA-1), Alexa Fluor® 647 anti-human CD15 (SSEA-1), PE/Cyanine5 anti-human CD15 (SSEA-1), PerCP anti-human CD15 (SSEA-1), PerCP/Cyanine5.5 anti-human CD15 (SSEA-1), Pacific Blue™ anti-human CD15 (SSEA-1), Alexa Fluor® 700 anti-human CD15 (SSEA-1), Brilliant Violet 510™ anti-human CD15 (SSEA-1), PE/Cyanine7 anti-human CD15 (SSEA-1), Brilliant Violet 605™ anti-human CD15 (SSEA-1), Brilliant Violet 650™ anti-human CD15 (SSEA-1), Purified anti-human CD15 (SSEA-1) (Maxpar® Ready), PE/Dazzle™ 594 anti-human CD15 (SSEA-1), Brilliant Violet 421™ anti-human CD15 (SSEA-1), APC/Fire™ 750 anti-human CD15 (SSEA-1), TotalSeq™-A0392 anti-human CD15 (SSEA-1), Brilliant Violet 711™ anti-human CD15 (SSEA-1), APC/Cyanine7 anti-human CD15 (SSEA-1), TotalSeq™-B0392 anti-human CD15 (SSEA-1), TotalSeq™-C0392 anti-human CD15 (SSEA-1), APC/Fire™ 810 anti-human CD15 (SSEA-1) Antibody, TotalSeq™-D0392 anti-human CD15 (SSEA-1), GMP FITC anti-human CD15 (SSEA-1), Spark YG™ 593 anti-human CD15 (SSEA-1), GMP PE anti-human CD15 (SSEA-1), GMP Pacific Blue™ anti-human CD15 (SSEA-1)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587