

## LEGENDplex™ Human Apolipoprotein (Apo) H Capture Bead B6, 13X

|                          |   |
|--------------------------|---|
| <b>Catalog# / Size</b>   | 740465 / 100 tests (270 µL)   |
| <b>Regulatory Status</b> | RUO   |
| <b>Other Names</b>       | Apolipoprotein H, Activated Protein C-Binding Protein, Anticardiolipin Cofactor, Beta-2-Glycoprotein I, APC Inhibitor   |
| <b>Description</b>       | <p>LEGENDplex™ Human Apolipoprotein (Apo) Panel Capture Beads are intended for use with the following reagents:</p> <p>740456 (LEGENDplex™ Human Apolipoprotein (Apo) Panel Detection Antibodies)<br/>740455 (LEGENDplex™ Human Apolipoprotein (Apo) Panel Standard)<br/>740375 (LEGENDplex™ Buffer Set D)</p> <p>The assembled Human Apolipoprotein (Apo) Mix and Match Panel is a multiplex bead-based assay panel, using fluorescence-encoded beads suitable for use on various flow cytometers. This allows simultaneous quantification of up to 11 human Apo proteins, including Apo AI, All, B100, CII, CIII, D, E, E4, H, J and M. This assay panel provides broader dynamic ranges than traditional ELISA methods. These targets have been validated for use with serum, plasma and cell culture supernatant samples.</p> |

### Product Details

---

|                            |  |
|----------------------------|--|
| <b>Verified Reactivity</b> | Human  |
| <b>Application</b>         | <a href="#">Multiplex</a>  |
|                            | Learn more about LEGENDplex™ at <a href="http://biolegend.com/legendplex">biolegend.com/legendplex</a> |
|                            | Download the <a href="#">LEGENDplex™ software here</a> .   |

### Antigen Details

---

|                         |                        |
|-------------------------|------------------------|
| <b>Biology Area</b>     | Immunology             |
| <b>Molecular Family</b> | Enzymes and Regulators |
| <b>Gene ID</b>          | <a href="#">350</a>    |

### Related Protocols

---

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587