

## LEGENDplex™ Human OPN Capture Bead A10, 13X

<b>Catalog# / Size</b>	740560 / 100 tests (270 µL)
<b>Regulatory Status</b>	RUO
<b>Other Names</b>	Osteopontin, SPP1, Secreted Phosphoprotein 1, Early T-Lymphocyte Activation 1, Urinary Stone Protein, Nephropontin, Uropontin, BNSP, Bone Sialoprotein I, ETA-1
<b>Description</b>	<p>The Vascular Inflammation Panel 1 Mix and Match Capture Beads offer flexible customization within the panel (13-plex) with the following items listed below:</p> <ul style="list-style-type: none"><li>740555 LEGENDplex™ Human Myoglobin Capture Bead A4, 13X</li><li>740556 LEGENDplex™ Human MRP8/14 Capture Bead A5, 13X</li><li>740557 LEGENDplex™ Human NGAL Capture Bead A6, 13X</li><li>740558 LEGENDplex™ Human CRP Capture Bead A7, 13X</li><li>740559 LEGENDplex™ Human MMP-2 Capture Bead A8, 13X</li><li>740560 LEGENDplex™ Human OPN Capture Bead A10, 13X</li><li>740561 LEGENDplex™ Human MPO Capture Bead B2, 13X</li><li>740562 LEGENDplex™ Human SAA Capture Bead B3, 13X</li><li>740476 LEGENDplex™ Human IGFBP-4 Capture Bead B4, 13X</li><li>740563 LEGENDplex™ Human ICAM-1 Capture Bead B5, 13X</li><li>740564 LEGENDplex™ Human VCAM-1 Capture Bead B6, 13X</li><li>740565 LEGENDplex™ Human MMP-9 Capture Bead B7, 13X</li><li>740566 LEGENDplex™ Human Cystatin C Capture Bead B9, 13X</li></ul>

Selected items will be used together with:

740554 LEGENDplex™ Human Vascular Inflammation Panel 1 Standard  
740553 LEGENDplex™ Human Vascular Inflammation Panel 1 Detection Antibodies  
740374 LEGENDplex™ Buffer Set C  
740377 or 740379 Filter Plate or V-bottom Plate

The assembled Human Vascular Inflammation Mix and Match Panel 1 is a multiplex bead-based assay kit, using fluorescence–encoded beads suitable for use on various flow cytometers. This allows simultaneous quantification of up to 13 human proteins: Myoglobin, MRP8/14, NGAL, CRP, MMP-2, OPN, MPO, SAA, IGFBP-4, ICAM-1, VCAM-1, MMP-9, and Cystatin C. This assay panel provides higher detection sensitivities and broader dynamic ranges than traditional ELISA methods. The panel has been validated for use with serum, plasma, tissue culture supernatant, and urine samples.

### Product Details

---

<b>Verified Reactivity</b>	Human
<b>Application</b>	<a href="#">Multiplex</a>
	Learn more about LEGENDplex™ at <a href="http://biolegend.com/legendplex">biolegend.com/legendplex</a>
	Download the <a href="#">LEGENDplex™ software here</a> .

### Antigen Details

---

<b>Molecular Family</b>	Enzymes and Regulators
<b>Gene ID</b>	<a href="#">6696</a>

### Related Protocols

---

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587