

## TotalSeq™-A0095 Rat IgG2b, κ Isotype Ctrl Antibody

<b>Catalog# / Size</b>	400673 / 10 µg
<b>Clone</b>	RTK4530
<b>Regulatory Status</b>	RUO
<b>Isotype</b>	Rat IgG2b, κ
<b>Barcode Sequence</b>	GATTCTTGACGACCT
<b>Description</b>	The isotype of RTK4530 immunoglobulin is rat IgG2b, κ. This antibody was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse, rat and human tissues.

### Product Details

<b>Antibody Type</b>	Monoclonal
<b>Host Species</b>	Rat
<b>Immunogen</b>	Trinitrophenol + KLH
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 1 mM EDTA.
<b>Preparation</b>	The antibody was purified by chromatography and conjugated with TotalSeq™-A oligomer under optimal conditions.
<b>Concentration</b>	0.5 mg/ml
<b>Storage &amp; Handling</b>	The antibody solution should be stored undiluted between 2°C and 8°C. <b>Do not freeze.</b>
<b>Recommended Usage</b>	<p>Each lot of this antibody is quality control tested by <a href="#">immunofluorescent staining with flow cytometric analysis</a> and the oligomer sequence is confirmed by sequencing. TotalSeq™-A antibodies are compatible with 10x Genomics Single Cell Gene Expression <a href="#">Solutions</a>.</p> <p>To maximize performance, it is strongly recommended that the reagent be titrated for each application, and that you centrifuge the antibody dilution before adding to the cells at 14,000xg at 2 - 8°C for 10 minutes. Carefully pipette out the liquid avoiding the bottom of the tube and add to the cell suspension. For Proteogenomics analysis, the suggested starting amount of this reagent for titration is ≤ 1.0 µg per million cells in 100 µL volume. Refer to the corresponding TotalSeq™ protocol for specific staining instructions.</p> <p>Buyer is solely responsible for determining whether Buyer has all intellectual property rights that are necessary for Buyer's intended uses of the BioLegend TotalSeq™ products. For example, for any technology platform Buyer uses with TotalSeq™, it is Buyer's sole responsibility to determine whether it has all necessary third party intellectual property rights to use that platform and TotalSeq™ with that platform.</p>
<b>Application Notes</b>	The RTK4530 immunoglobulin is useful as an isotype-matched control (for the relevant formats) for Western blotting, immunoprecipitation, immunohistochemistry, functional assay, and immunofluorescence microscopy. The Ultra-LEAF™ purified antibody (Endotoxin < 0.01 EU/µg, Azide-Free, 0.2 µm filtered) is recommended for functional assays (Cat. No. 400643, 400644, 400671, 400672, 400675, and 400676) as negative control.
<b>Additional Product Notes</b>	<p>TotalSeq™ reagents are designed to profile protein levels at a single cell level following an optimized protocol similar to the CITE-seq workflow. A compatible single cell device (e.g. <a href="#">10x Genomics Chromium System and Reagents</a>) and sequencer (e.g. Illumina analyzers) are required. Please contact <a href="#">technical support</a> for more information, or visit <a href="#">biolegend.com/totalseq</a>.</p> <p>The barcode flanking sequences are CCTTGGCACCCGAGAATTCCA (PCR handle), and BAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA*A*A (capture sequence). B represents either C, G, or T, and * indicates a phosphorothioated bond, to prevent nuclease degradation.</p>

View more applications data for this product in our [Scientific Poster Library](#).

### Application References

**(PubMed link indicates BioLegend citation)**

1. Cervantes-Barragan L, *et al.* 2007. *Blood* 109:1131.
2. Zeiser R, *et al.* 2007. *Blood* 109:2225.
3. Sasaki K, *et al.* 2008. *J. Immunol.* 181:104. [PubMed](#)
4. Duan J, *et al.* 2008. *P. Natl. Acad. Sci. USA* 105:5183. [PubMed](#)
5. Yi H, *et al.* 2009. *Blood* 113:5819. [PubMed](#)
6. Schafeer JS, *et al.* 2010. *J. Leukocyte Biol.* 87:301. [PubMed](#)
7. Lei GS, *et al.* 2015. *Infect Immun.* 83:572. [PubMed](#)
8. Richards J, *et al.* 2015. *Mol Cell Cardiol.* 79:21. [PubMed](#)

**Product Citations**

1. Hao Y, *et al.* 2021. *Cell.* 184:3573. [PubMed](#)

**Antigen Details**

---

<b>Structure</b>	The antibody was purified by affinity chromatography and conjugated with APC/Fire™ 810 under optimal conditions.
<b>Gene ID</b>	NA

**Related Protocols**

---

[TotalSeq™-A Antibodies and Cell Hashing with 10x Single Cell 3' Reagent Kit v3 3.1 Protocol](#)

**Other Formats**

---

APC Rat IgG2b, κ Isotype Ctrl, Biotin Rat IgG2b, κ Isotype Ctrl, FITC Rat IgG2b, κ Isotype Ctrl, PE Rat IgG2b, κ Isotype Ctrl, PE/Cyanine5 Rat IgG2b, κ Isotype Ctrl, Purified Rat IgG2b, κ Isotype Ctrl, PE/Cyanine7 Rat IgG2b, κ Isotype Ctrl, APC/Cyanine7 Rat IgG2b, κ Isotype Ctrl, Alexa Fluor® 647 Rat IgG2b, κ Isotype Ctrl, Alexa Fluor® 488 Rat IgG2b, κ Isotype Ctrl, Pacific Blue™ Rat IgG2b, κ Isotype Ctrl, Alexa Fluor® 700 Rat IgG2b, κ Isotype Ctrl, PerCP Rat IgG2b, κ Isotype Ctrl, PerCP/Cyanine5.5 Rat IgG2b, κ Isotype Ctrl, Brilliant Violet 421™ Rat IgG2b, κ Isotype Ctrl, Brilliant Violet 570™ Rat IgG2b, κ Isotype Ctrl, Ultra-LEAF™ Purified Rat IgG2b, κ Isotype Ctrl, Brilliant Violet 510™ Rat IgG2b, κ Isotype Ctrl, Brilliant Violet 605™ Rat IgG2b, κ Isotype Ctrl, Brilliant Violet 785™ Rat IgG2b, κ Isotype Ctrl, Brilliant Violet 650™ Rat IgG2b, κ Isotype Ctrl, Brilliant Violet 711™ Rat IgG2b, κ Isotype Ctrl, PE/Dazzle™ 594 Rat IgG2b, κ Isotype Ctrl, Alexa Fluor® 594 Rat IgG2b, κ Isotype Ctrl, GolnVivo™ Purified Rat IgG2b, κ Isotype Ctrl, APC/Fire™ 750 Rat IgG2b, κ Isotype Ctrl, TotalSeq™-A0095 Rat IgG2b, κ Isotype Ctrl, TotalSeq™-C0095 Rat IgG2b, κ Isotype Ctrl, KIRAVIA Blue 520™ Rat IgG2b, κ Isotype Ctrl, Spark Blue™ 550 Rat IgG2b, κ Isotype Ctrl, PE/Fire™ 640 Rat IgG2b, κ Isotype Ctrl, APC/Fire™ 810 Rat IgG2b, κ Isotype Ctrl, PE/Fire™ 700 Rat IgG2b, κ Isotype Ctrl, TotalSeq™-B0095 Rat IgG2b, κ Isotype Ctrl

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587