

## LEGENDplex™ Human FGF-basic Capture Bead A7, 13X

<b>Catalog# / Size</b>	740184 / 100 tests (270 µL)
<b>Regulatory Status</b>	RUO
<b>Description</b>	LEGENDplex™ Human Growth Factor Mix and Match Panel Capture Beads are intended for use with the following reagents:  740195 (LEGENDplex™ Human Growth Factor Detection Antibodies) 740194 (LEGENDplex™ Human Growth Factor Standard) 741148 (LEGENDplex™ Buffer Set A5)

### Product Details

---

<b>Verified Reactivity</b>	Human
<b>Application</b>	<a href="#">Multiplex</a>  Learn more about LEGENDplex™ at <a href="http://biolegend.com/legendplex">biolegend.com/legendplex</a>  Download the <a href="#">LEGENDplex™ software here</a> .

### Antigen Details

---

<b>Biology Area</b>	Cell Biology, Neuroscience, Synaptic Biology
<b>Molecular Family</b>	Cytokines/Chemokines, Growth Factors
<b>Gene ID</b>	<a href="#">2247</a>

### Related Protocols

---

[LEGENDplex™ Laboratory Protocol - Video](#)

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587