

LEGENDplex™ Human sCD27 Capture Bead A7, 13X

Catalog# / Size	740873 / 100 tests (270 µL)
Regulatory Status	RUO
Other Names	T14, S152, Tp55, TNFRSF7, S152, LPFS2, CD27, sCD27, CD27 molecule
Description	<p>The LEGENDplex™ HU Immune Checkpoint Panel 1 Mix and Match Beads are individual capture beads for 12-targets including: sCD25 (IL-2Ra), 4-1BB, sCD27, B7.2 (CD86), Free Active TGF-b1, CTLA-4, PD-L1, PD-L2, PD-1, Tim-3, LAG-3, and Galectin-9. The beads are intended for use with each other and the following in Mix and Match Assays:</p> <p>740868 LEGENDplex™ HU Immune Checkpoint Panel 1 Detection Antibodies 740869 LEGENDplex™ HU Immune Checkpoint Panel 1 Standard 740623 LEGENDplex™ Buffer Set M 740377 LEGENDplex™ Filter Plate OR 740379 LEGENDplex™ V-bottom Plate</p>

Product Details

Verified Reactivity	Human
Application	Multiplex
	Learn more about LEGENDplex™ at biolegend.com/legendplex
	Download the LEGENDplex™ software here .

Antigen Details

Biology Area	Immunology
Molecular Family	CD Molecules
Gene ID	939

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 www.biolegend.com
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587