

## TotalSeq™-C0236 Rat IgG1, κ Isotype Ctrl Antibody

<b>Catalog# / Size</b>	400467 / 10 µg
<b>Clone</b>	RTK2071
<b>Regulatory Status</b>	RUO
<b>Isotype</b>	Rat IgG1, κ
<b>Barcode Sequence</b>	ATCAGATGCCCTCAT
<b>Description</b>	The isotype of RTK2071 immunoglobulin is rat IgG1, κ. This antibody was chosen as an isotype control after screening on a variety of resting, activated, live, and fixed mouse, rat and human tissues.

### Product Details

<b>Antibody Type</b>	Monoclonal
<b>Host Species</b>	Rat
<b>Immunogen</b>	Trinitrophenol + KLH
<b>Formulation</b>	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 1 mM EDTA
<b>Preparation</b>	The antibody was purified by chromatography and conjugated with TotalSeq™-C oligomer under optimal conditions.
<b>Concentration</b>	0.5 mg/mL
<b>Storage &amp; Handling</b>	The antibody solution should be stored undiluted between 2°C and 8°C. <b>Do not freeze.</b>
<b>Application</b>	<a href="#">PG - Quality tested</a>
<b>Recommended Usage</b>	<p>Each lot of this antibody is quality control tested by <a href="#">immunofluorescent staining with flow cytometric analysis</a> and the oligomer sequence is confirmed by sequencing. TotalSeq™-C antibodies are compatible with 10x Genomics Chromium Single Cell Immune Profiling <a href="#">Solution</a>.</p> <p>To maximize performance, it is strongly recommended that the reagent be titrated for each application, and that you centrifuge the antibody dilution before adding to the cells at 14,000xg at 2 - 8°C for 10 minutes. Carefully pipette out the liquid avoiding the bottom of the tube and add to the cell suspension. For Proteogenomics analysis, the suggested starting amount of this reagent for titration is ≤ 1.0 µg per million cells in 100 µL volume. Refer to the corresponding TotalSeq™ protocol for specific staining instructions.</p> <p>Buyer is solely responsible for determining whether Buyer has all intellectual property rights that are necessary for Buyer's intended uses of the BioLegend TotalSeq™ products. For example, for any technology platform Buyer uses with TotalSeq™, it is Buyer's sole responsibility to determine whether it has all necessary third party intellectual property rights to use that platform and TotalSeq™ with that platform.</p>
<b>Application Notes</b>	The RTK2071 immunoglobulin is useful as an isotype-matched control (for the relevant formats) for Western blotting, immunoprecipitation, immunohistochemistry, functional assay, and immunofluorescence microscopy. The LEAF™ purified antibody (Endotoxin < 0.1 EU/µg, Azide-Free, 0.2 µm filtered) is recommended for functional assays (Cat. No. 400414) as negative control <sup>1</sup> . For <i>in vivo</i> studies or highly sensitive assays, we recommend Ultra-LEAF™ purified antibody (Cat. No. 400432) with a lower endotoxin limit than standard LEAF™ purified antibodies (Endotoxin < 0.01 EU/µg).
<b>Additional Product Notes</b>	<p>TotalSeq™ reagents are designed to profile protein levels at a single cell level following an optimized protocol similar to the CITE-seq workflow. A compatible single cell device (e.g. <a href="#">10x Genomics Chromium System and Reagents</a>) and sequencer (e.g. Illumina analyzers) are required. Please contact <a href="#">technical support</a> for more information, or visit <a href="http://biolegend.com/totalseq">biolegend.com/totalseq</a>.</p> <p>The barcode flanking sequences are CGGAGATGTGTATAAGAGACAGNNNNNNNNNN (PCR handle), and NNNNNNNNCCCATATAAGA*A*A (capture sequence). N represents either randomly selected A, C, G, or T, and * indicates a phosphorothioated bond, to prevent nuclease</p>

degradation.

View more applications data for this product in our [Scientific Poster Library](#).

## Application References

(PubMed link indicates BioLegend citation)

1. Riemann M, *et al.* 2005. *J. Immunol.* 175:3560. [PubMed](#)
2. Wondimu Z, *et al.* 2010. *Am. J. Pathol.* 177:2334. [PubMed](#)

## Antigen Details

---

Gene ID                      NA

## Related Protocols

---

[TotalSeq™-B or -C with 10x Feature Barcoding Technology](#)

## Other Formats

---

APC Rat IgG1, κ Isotype Ctrl, Biotin Rat IgG1, κ Isotype Ctrl, FITC Rat IgG1, κ Isotype Ctrl, PE Rat IgG1, κ Isotype Ctrl, PE/Cyanine5 Rat IgG1, κ Isotype Ctrl, Purified Rat IgG1, κ Isotype Ctrl, PE/Cyanine7 Rat IgG1, κ Isotype Ctrl, Alexa Fluor® 647 Rat IgG1, κ Isotype Ctrl, Alexa Fluor® 488 Rat IgG1, κ Isotype Ctrl, Pacific Blue™ Rat IgG1, κ Isotype Ctrl, Alexa Fluor® 700 Rat IgG1, κ Isotype Ctrl, APC/Cyanine7 Rat IgG1, κ Isotype Ctrl, PerCP Rat IgG1, κ Isotype Ctrl, PerCP/Cyanine5.5 Rat IgG1, κ Isotype Ctrl, Brilliant Violet 421™ Rat IgG1, κ Isotype Ctrl, Ultra-LEAF™ Purified Rat IgG1, κ Isotype Ctrl, Brilliant Violet 605™ Rat IgG1, κ Isotype Ctrl, Brilliant Violet 510™ Rat IgG1, κ Isotype Ctrl, Brilliant Violet 650™ Rat IgG1, κ Isotype Ctrl, Brilliant Violet 711™ Rat IgG1, κ Isotype Ctrl, Brilliant Violet 785™ Rat IgG1, κ Isotype Ctrl, PE/Dazzle™ 594 Rat IgG1, κ Isotype Ctrl, Alexa Fluor® 594 Rat IgG1, κ Isotype Ctrl, GolnVivo™ Purified Rat IgG1, κ Isotype Ctrl, APC/Fire™ 750 Rat IgG1, κ Isotype Ctrl, TotalSeq™-A0236 Rat IgG1, κ Isotype Ctrl, KIRAVIA Blue 520™ Rat IgG1, κ Isotype Ctrl, TotalSeq™-B0236 Rat IgG1, κ Isotype Ctrl, TotalSeq™-C0236 Rat IgG1, κ Isotype Ctrl, Spark YG™ 581 Rat IgG1, κ Isotype Ctrl, Spark YG™ 593 Rat IgG1, κ Isotype Ctrl

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

\*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, [www.biolegend.com/ordering#license](http://www.biolegend.com/ordering#license)). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.

BioLegend Inc., 8999 BioLegend Way, San Diego, CA 92121 [www.biolegend.com](http://www.biolegend.com)  
Toll-Free Phone: 1-877-Bio-Legend (246-5343) Phone: (858) 768-5800 Fax: (877) 455-9587