

LEGENDplex™ Human IL-4 Capture Bead A4, 13X

Catalog# / Size	740933 / 100 tests (270 µL)
Regulatory Status	RUO
Other Names	IL4, BSF1, IL-4, BCGF1, BSF-1, BCGF-1, Interleukin 4
Description	<p>The LEGENDplex™ Human Essential Immune Response Panel Mix and Match Beads are individual capture beads for 13 targets including IL-4, IL-2, CXCL10 (IP-10), IL-1β, TNF-α, CCL2 (MCP-1), IL-17A, IL-6, IL-10, IFN-γ, IL-12p70, CXCL8 (IL-8), Free Active TGF-β1. The beads are recommended for use with each other and the following in Mix and Match assays:</p> <p>Human Essential Immune Response Panel</p> <p>740932 LEGENDplex™ Human Essential Immune Response Panel Standard 740931 LEGENDplex™ Human Essential Immune Response Panel Detection Antibodies 740368 LEGENDplex™ Buffer Set A 740377 LEGENDplex™ Filter Plate OR 740379 LEGENDplex™ V-Bottom Plate</p>

Product Details

Verified Reactivity	Human
Application	Multiplex Learn more about LEGENDplex™ at biolegend.com/legendplex Download the LEGENDplex™ software here .
Product Citations	1. Schiferle EB, <i>et al.</i> 2021. Sci Adv. 7: . PubMed

Antigen Details

Biology Area	Immunology
Molecular Family	Cytokines/Chemokines
Gene ID	3565

For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.

*These products may be covered by one or more Limited Use Label Licenses (see the BioLegend Catalog or our website, www.biolegend.com/ordering#license). BioLegend products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products, reverse engineer functionally similar materials, or to provide a service to third parties without written approval of BioLegend. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.